

Disinfectants for SonoSite Products

Table of Contents

Introduction	1
Ultrasound Systems and Mini-Docks	2
Standard Transducers, I-Z	12
Accessories	20
Cleaning transducers	25

Introduction

SonoSite tests disinfectants and disinfectant devices for use with its systems, transducers, and accessories. Except for TEE and surgical transducers, results are in the tables below. See the ultrasound system user guide for cleaning and disinfecting instructions and for additional warnings and cautions.

For cleaners and disinfectants compatible with TEE and surgical transducers, see the *TEE Care Card* and *Surgical Transducer Guide*.

WARNING: Disinfectants listed are recommended by SonoSite for safe use with SonoSite products.
Refer to instructions on the disinfectant label for guidance on disinfection effectiveness and appropriate clinical use.

This document does not have the following regulatory information for disinfectants:

- EPA registration
- FDA 510(k) clearance (liquid sterilant, high level disinfectant)
- CE approval

Before using a disinfectant, confirm that its regulatory status is appropriate for your jurisdiction and use.

Verify expiration dates on chemicals.

When disposing of chemicals, follow manufacturer recommendations and EPA regulations.

Ultrasound Systems and Mini-Docks

You can disinfect the exterior surface of the ultrasound system using a recommended disinfectant.

Table 1: Disinfectant Compatibility with Sonosite Ultrasound Systems and Mini-Docks

Disinfectant	Edge™ (system & mini-dock)	NanoMaxx®	S Series™	M-Turbo® M-OB/GYN Office™ (system & mini-dock)	MicroMaxx® (system & mini-dock)	TITAN® (system & mini-dock)	SonoSite® 180PLUS™	iLook®15 iLook®25
AbcoCide 14								✓
Airkem A-33®	✓	✓	✓	✓	✓	✓	✓	✓
Alcohol, Ethanol	✓	✓	✓					
Alkacide							✓	
Alkazyme							✓	
Anios Wipes	✓	✓		✓	✓			✓
Anioxyde 1000	✓				✓		✓	
Aniosyme DD1 (5%)	✓	✓	✓					
Aquatabs (1000 PPM)							✓	
Ascend							✓	✓
Asepti-HB	✓	✓	✓	✓				
Asepti-Steryl®	✓							
Asepti-Wipes						✓		
Asepti-Wipe™ II	✓		✓					
Bacillocid rasant							✓	✓

Table 1: Disinfectant Compatibility with Sonosite Ultrasound Systems (continued) and Mini-Docks

Disinfectant	Edge™ (system & mini-dock)	NanoMaxx®	S Series™	M-Turbo® M-OB/GYN Office™ (system & mini-dock)	MicroMaxx® (system & mini-dock)	TITAN® (system & mini-dock)	SonoSite® 180PLUS™	iLook®15 iLook®25
Banicide								✓
Bacoban WB			✓				✓	
Baxedin Pre-Op	✓	✓	✓		✓	✓		
BioText Micrylium							✓	
CaviWipes™	✓				✓	✓	✓	
Cidalkan®	✓			✓	✓	✓		
Cidalkan® Lingettes	✓	✓			✓	✓	✓	
Cidex®	✓	✓	✓	✓	✓	✓		
Cidex® 14							✓	✓
Cidex® OPA	✓	✓	✓	✓		✓		✓
Cidex® Plus	✓	✓	✓	✓	✓	✓	✓	✓
Cleanisept Wipe			✓					✓
Clinell Universal Sanitising Wipes	✓	✓			✓	✓		
Clorox Germicidal Wipes						✓		
Control III	✓	✓	✓		✓	✓		
Coverage Germicidal Surface Wipes	✓	✓						
Coverage Plus Germicidal Surface Wipes	✓	✓				✓		

Table 1: Disinfectant Compatibility with Sonosite Ultrasound Systems (continued)and Mini-Docks

Disinfectant	Edge™ (system & mini-dock)	NanoMaxx®	S Series™	M-Turbo® M-OB/GYN Office™ (system & mini-dock)	MicroMaxx® (system & mini-dock)	TITAN® (system & mini-dock)	SonoSite® 180PLUS™	iLook®15 iLook®25
Coverage™ Spray	✓							
DentaSept							✓	
Discide Ultra Wipes	✓		✓					
DisOPA								✓
Dispatch Wipes					✓			
Dynacide PA	✓	✓	✓	✓	✓	✓	✓	✓
Echo Clean Lingettes	✓	✓	✓	✓				
End-Bac® II	✓				✓	✓		
Endosporine	✓	✓		✓	✓			✓
EndoZime AW® Plus							✓	✓
Envirocide®	✓		✓					
Enzol							✓	✓
Expose™ II							✓	
Gluco Chlor Wipes					✓			
Hexanios							✓	✓
Hi Tor Plus							✓	
Hibiclens							✓	✓

Table 1: Disinfectant Compatibility with Sonosite Ultrasound Systems (continued)and Mini-Docks

Disinfectant	Edge™ (system & mini-dock)	NanoMaxx®	S Series™	M-Turbo® M-OB/GYN Office™ (system & mini-dock)	MicroMaxx® (system & mini-dock)	TITAN® (system & mini-dock)	SonoSite® 180PLUS™	iLook®15 iLook®25
Hydrogen Peroxide (3%)	✓	✓	✓	✓				
Incidin Plus (1%)	✓	✓	✓		✓	✓	✓	
Incidin Plus (3%)	✓	✓			✓		✓	
Korsolex® Extra							✓	✓
MadaCide 1	✓		✓					
Medzyme	✓	✓	✓	✓	✓	✓		
MetriCide 14							✓	✓
MetriCide® 28							✓	✓
MetriCide OPA Plus	✓	✓		✓		✓		✓
MetriZyme							✓	✓
Mikrobac forte							✓	✓
Mikrobac forte 2%	✓							
Mikrobac Tissues	✓							
Nu-Cidex			✓		✓	✓		
Nuclean							✓	✓
PerCept RTU Wipes™		✓	✓	✓	✓	✓		
Rely+On™ PeraSafe™		✓				✓		✓

Table 1: Disinfectant Compatibility with Sonosite Ultrasound Systems (continued) and Mini-Docks

Disinfectant	Edge™ (system & mini-dock)	NanoMaxx®	S Series™	M-Turbo® M-OB/GYN Office™ (system & mini-dock)	MicroMaxx® (system & mini-dock)	TITAN® (system & mini-dock)	SonoSite® 180PLUS™	iLook®15 iLook®25
Salvanios pH 7					✓		✓	✓
Sani-Cloth® AF	✓							
Sani-Cloth® HB	✓	✓	✓	✓				
Sani-Cloth® Bleach Wipes		✓			✓			
Sani-Cloth® Plus	✓	✓	✓	✓	✓			✓
Sekusept® aktiv						✓		
Simple Chlor		✓	✓	✓				
Sporicidin®	✓						✓	
Steranios 2%	✓	✓	✓	✓	✓	✓		✓
Steranios 20%	✓	✓	✓		✓	✓		✓
Super Sani-Cloth®	✓							
T-Spray					✓	✓		
T-Spray II™	✓	✓	✓	✓	✓	✓	✓	✓
TASK 105							✓	✓
Trigene Advance Wipes	✓	✓			✓	✓		✓
Tristel Solo	✓	✓	✓		✓	✓	✓	✓
Tristel Wipes			✓	✓	✓	✓		

Table 1: Disinfectant Compatibility with Sonosite Ultrasound Systems (continued) and Mini-Docks

Disinfectant	Edge™ (system & mini-dock)	NanoMaxx®	S Series™	M-Turbo® M-OB/GYN Office™ (system & mini-dock)	MicroMaxx® (system & mini-dock)	TITAN® (system & mini-dock)	SonoSite® 180PLUS™	iLook®15 iLook®25
Virex® 256		✓						
Virex® II 256	✓	✓					✓	
Virex® TB	✓	✓						
Virox® 5 Wipes	✓	✓	✓		✓	✓	✓	
Virufen							✓	✓
Wavicide -01								✓
Wavicide -06							✓	
Wet Task Wipes	✓	✓	✓	✓	✓	✓		
Wet Wipe Disinfection					✓	✓	✓	✓
Wex-Cide	✓	✓	✓	✓		✓	✓	✓
Wipes Plus Disinfecting Wipes			✓		✓	✓		
Wood Wyant Ultra Wipes	✓		✓					
✓ - Acceptable for use								

Standard Transducers, A-H

You can disinfect the transducer, scanning surface, and cable using a recommended disinfectant or device. Clean the transducer prior to disinfecting to remove any particulate matter or bodily fluids. See [“To clean a transducer”](#) on page 25.

Table 2: Disinfectant-Device Compatibility with SonoSite Standard Transducers (A-H)

Device	C11 C11e C11x	C11n	C11 Cast Urethane ¹	C15 C15e	C15 Cast Urethane ²	C60 C60e C60x	C60n	D2 D2x	HFL38 HFL38x	HFL50x
Nanosonics Trophon EPR*	✓	✓		✓		✓	✓		✓	✓

✓ - Acceptable for use

1- Part number: P01960-xx, P02645-02, P02645-03, or P03361-10

2 - Part number: P00798-03

* For information and guidance on disinfection effectiveness and appropriate clinical use, see the manufacturer's user manual.

Table 3: Disinfectant Compatibility with SonoSite Standard Transducers (A-H)

Disinfectant	C11 C11e C11x	C11n	C11 Cast Urethane ¹	C15 C15e	C15 Cast Urethane ²	C60 C60e C60x	C60n	D2 D2x	HFL38 HFL38x	HFL50x
Airkem A-33®	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Alcohol, Ethanol		✓					✓			✓
Anios Wipes	✓		✓	✓	✓	✓		✓	✓	✓
Aniosyme DD1 (5%)	✓	✓	✓	✓	✓	✓	✓		✓	✓
Anioxyde 1000								✓		
Asepti-HB		✓					✓			✓
Asepti-Wipe™ II	✓		✓	✓	✓	✓			✓	
Bacoban	✓		✓	✓	✓	✓		✓	✓	
Bacoban WB								✓		

Table 3: Disinfectant Compatibility with SonoSite Standard Transducers (A-H) (continued)

Disinfectant	C11 C11e C11x	C11n	C11 Cast Urethane 1	C15 C15e	C15 Cast Urethane 2	C60 C60e C60x	C60n	D2 D2x	HFL38 HFL38x	HFL50x
Baxedin Pre-Op	✓	✓	✓	✓	✓	✓	✓		✓	✓
CaviWipes™	✓	✓	✓	✓	✓	✓	✓		✓	✓
Cidalkan®	✓		✓	✓	✓	✓			✓	
Cidalkan® Lingettes	✓		✓	✓	✓	✓			✓	
Cidex®	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Cidex® OPA	✓	✓ ³	✓	✓	✓	✓	✓ ³	✓	✓	✓
Cidex® Plus	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Cleanisept Wipes	✓		✓	✓	✓	✓		✓	✓	
Clinell Universal Sanitising Wipes	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Clorox® Wipes					✓					
Clorox® Germicidal Wipes	✓	✓	✓	✓	✓	✓	✓		✓	✓
Control III	✓		✓	✓	✓	✓				
Coverage™ Germicidal Surface Wipes	✓	✓	✓	✓	✓	✓	✓		✓	✓
Coverage™ Plus Germicidal Surface Wipes	✓	✓	✓	✓	✓	✓	✓		✓	✓
Coverage™ Spray		✓					✓			
DisOPA								✓		
Dispatch	✓		✓		✓					
Dispatch Wipes				✓		✓		✓	✓	
Dynacide PA	✓	✓	✓	✓	✓	✓	✓	✓	✓	

Table 3: Disinfectant Compatibility with SonoSite Standard Transducers (A-H) (continued)

Disinfectant	C11 C11e C11x	C11n	C11 Cast Urethane 1	C15 C15e	C15 Cast Urethane 2	C60 C60e C60x	C60n	D2 D2x	HFL38 HFL38x	HFL50x
Echo Clean Linettes	✓		✓	✓	✓	✓			✓	✓
End-Bac® II	✓		✓	✓	✓	✓			✓	✓
Endosporine	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Envirocide®		✓					✓			
Gluco Chlor Wipes	✓	✓	✓	✓	✓	✓	✓		✓	✓
Hydrogen Peroxide (3%)	✓	✓	✓	✓	✓	✓	✓	✓		✓
Incidin Plus (1%)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Incidin Plus (3%)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Korsolex® Extra	✓		✓	✓	✓	✓			✓	
Madacide 1	✓	✓	✓	✓	✓	✓	✓		✓	✓
Medzyme	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Metricide OPA Plus	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
MetriSponge	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Nu-Cidex								✓		
PerCept RTU Wipes™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rely+On™ PeraSafe™								✓		
Sani-Cloth® Bleach Wipes	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Sani-Cloth® HB	✓	✓	✓	✓	✓	✓	✓			✓
Sani-Cloth® Plus	✓		✓	✓	✓	✓				✓
Sekusept®	✓		✓	✓	✓	✓		✓	✓	

Table 3: Disinfectant Compatibility with SonoSite Standard Transducers (A-H) (continued)

Disinfectant	C11 C11e C11x	C11n	C11 Cast Urethane 1	C15 C15e	C15 Cast Urethane 2	C60 C60e C60x	C60n	D2 D2x	HFL38 HFL38x	HFL50x
Staphene®	✓		✓	✓	✓	✓			✓	✓
Steranios (2%)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Steranios 20%								✓		
T-Spray™		✓		✓		✓	✓			✓
T-Spray II™	✓	✓ ³	✓	✓	✓	✓	✓ ³	✓	✓	✓
Trigene Advance Wipes	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Tristel Duo								✓		
Tristel Solo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Tristel Wipes								✓		
Virex® 256		✓					✓			✓
Virex® TB		✓					✓			✓
Virox® 5 Wipes	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Virufen								✓		
Wet Task Wipes	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Wet Wipe Disinfection	✓		✓	✓	✓	✓		✓	✓	
Wex-Cide	✓		✓	✓	✓	✓				
Wipes Plus Disinfecting Wipes	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

✓ - Acceptable for use

1- Part number: P01960-xx, P02645-02, P02645-03, or P03361-10

2 - Part number: P00798-03

3 - Disinfectant has been tested by SonoSite for disinfection effectiveness on SonoSite transducers.

Standard Transducers, I-Z

You can disinfect the transducer, scanning surface, and cable using a recommended disinfectant or device. Clean the transducer prior to disinfecting to remove any particulate matter or bodily fluids. See “[To clean a transducer](#)” on page 25.

Table 4: Disinfectant-Device Compatibility with SonoSite Standard Transducers (I-Z)

Device	L25 L25e L25x	L25n	L38 L38e L38x	L38n L38xi	L52 L52e Vet Only	L52x Vet Only	L52n Vet Only	P10 P10x	P17	P21x	P21n
Nanosonics Trophon EPR*	✓	✓	✓	✓				✓	✓	✓	✓

✓ - Acceptable for use

* For information and guidance on disinfection effectiveness and appropriate clinical use, see the manufacturer’s user manual.

Table 5: Disinfectant Compatibility with SonoSite Standard Transducers (I-Z)

Disinfectant	L25 L25e L25x	L25n	L38 L38e L38x	L38n	L38xi	L52 L52e Vet Only	L52x Vet Only	L52n Vet Only	P10 P10x	P17	P21x	P21n
Airkem A-33®	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Alcohol, Ethanol		✓		✓	✓		✓	✓				✓
Anios Wipes	✓		✓			✓	✓		✓	✓	✓	
Aniosyme DD1 (5%)	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Anioxyde 1000						✓						
Asepti-HB		✓		✓	✓		✓	✓				✓
Asepti-Wipe™ II	✓		✓			✓			✓	✓	✓	
Bacoban	✓		✓			✓			✓	✓	✓	
Baxedin Pre-Op	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
CaviWipes™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

Table 5: Disinfectant Compatibility with SonoSite Standard Transducers (I-Z) (continued)

Disinfectant	L25 L25e L25x	L25n	L38 L38e L38x	L38n	L38xi	L52 L52e Vet Only	L52x Vet Only	L52n Vet Only	P10 P10x	P17	P21x	P21n
Cidalkan®	✓		✓			✓			✓	✓	✓	
Cidalkan® Lingettes	✓		✓			✓			✓	✓	✓	
Cidex®	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Cidex® OPA	✓	✓ ¹	✓	✓ ¹	✓	✓	✓	✓	✓	✓	✓	✓ ¹
Cidex® Plus	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Cleanisept Wipes	✓		✓			✓			✓	✓	✓	
Clinell Universal Sanitising Wipes	✓		✓	✓	✓			✓	✓	✓	✓	✓
Clorox® Germicidal Wipes	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Control III	✓		✓			✓			✓	✓	✓	
Coverage™ Germicidal Surface Wipes	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Coverage™ Plus Germicidal Surface Wipes	✓		✓	✓	✓		✓	✓	✓	✓	✓	✓
Coverage™ Spray				✓	✓			✓				✓
Dispatch	✓											
Dispatch Wipes			✓			✓			✓	✓	✓	
Dynacide PA	✓		✓	✓	✓	✓		✓	✓	✓	✓	✓
Echo Clean Lingettes	✓		✓			✓	✓		✓	✓	✓	
End-Bac® II	✓		✓			✓	✓		✓	✓	✓	
Endosporine	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓
Envirocide®				✓	✓			✓				✓

Table 5: Disinfectant Compatibility with SonoSite Standard Transducers (I-Z) (continued)

Disinfectant	L25 L25e L25x	L25n	L38 L38e L38x	L38n	L38xi	L52 L52e Vet Only	L52x Vet Only	L52n Vet Only	P10 P10x	P17	P21x	P21n
Gluco Chlor Wipes	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Hydrogen Peroxide (3%)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Incidin Plus (1%)	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Incidin Plus (3%)	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Korsolex® Extra	✓		✓			✓			✓	✓	✓	
Madacide 1	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Medzyme	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Metricide OPA Plus	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
MetriSponge	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
PerCept RTU Wipes™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rely+On™ PeraSafe™						✓	✓					
Sani-Cloth® Bleach Wipes	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Sani-Cloth® HB	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Sani-Cloth® Plus	✓		✓			✓	✓		✓	✓	✓	
Sekusept®	✓		✓			✓			✓	✓	✓	
Staphene®	✓		✓			✓	✓		✓	✓	✓	
Steranios (2%)	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Steranios (20%)						✓						
T-Spray™			✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
T-Spray II™	✓	✓ ¹	✓	✓ ¹	✓	✓	✓	✓	✓	✓	✓	✓ ¹

Table 5: Disinfectant Compatibility with SonoSite Standard Transducers (I-Z) (continued)

Disinfectant	L25 L25e L25x	L25n	L38 L38e L38x	L38n	L38xi	L52 L52e Vet Only	L52x Vet Only	L52n Vet Only	P10 P10x	P17	P21x	P21n
Trigene Advance Wipes	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Tristel Solo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Tristel Wipes						✓						
Virex® 256				✓	✓		✓	✓				✓
Virex® II 256		✓										
Virex® TB				✓	✓		✓	✓				✓
Virox® 5 Wipes	✓		✓	✓	✓	✓		✓	✓	✓	✓	✓
Virox® 5 RTU Wipes		✓										
Wet Task Wipes	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Wet Wipe Disinfection	✓		✓			✓			✓	✓	✓	
Wex-Cide	✓		✓						✓	✓	✓	
Wipes Plus Disinfecting Wipes	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

✓ - Acceptable for use

1 Disinfectant has been tested by SonoSite for disinfection effectiveness on SonoSite transducers.

Specialty Transducers

You can disinfect the transducer, scanning surface, and cable using a recommended disinfectant or device.

Clean the transducer prior to disinfecting to remove any particulate matter or bodily fluids. See [“To clean a transducer”](#) on page 25.

Table 6: Disinfectant-Device Compatibility with SonoSite Specialty Transducers

Device	C8 C8e C8x	ICT ICTx	ICTn	SLA SLAx	SLT
Nanosonics Trophon EPR*	✓	✓	✓	✓	

✓ - Acceptable for use

* For information and guidance on disinfection effectiveness and appropriate clinical use, see the manufacturer's user manual.

Table 7: Disinfectant Compatibility with SonoSite Specialty Transducers

Disinfectant	C8 C8e C8x	ICT ICTx	ICTn	SLA SLAx	SLT
Airkem A-33®	✓	✓	✓	✓	✓
Alcohol, Ethanol			✓		
Anios Wipes	✓	✓		✓	✓
Aniosyme DD1 (5%)	✓	✓		✓	✓
Anioxyde 1000					✓
Asepti-HB			✓		
Asepti-Wipe™ II	✓	✓		✓	✓
Bacoban	✓	✓		✓	✓
Baxedin Pre-Op	✓	✓	✓	✓	✓
CaviWipes™	✓	✓	✓	✓	✓
Cidalkan®	✓	✓		✓	✓

Table 7: Disinfectant Compatibility with SonoSite Specialty Transducers (continued)

Disinfectant	C8 C8e C8x	ICT ICTx	ICTn	SLA SLAx	SLT
Cidalkan® Lingettes	✓	✓		✓	✓
Cidex®	✓	✓	✓	✓	✓
Cidex® OPA	✓	✓	✓	✓	✓
Cidex® Plus	✓	✓	✓	✓	✓
Cleanisept Wipe	✓	✓		✓	✓
Clinell Universal Sanitising Wipes		✓	✓	✓	
Clorox® Germicidal Wipes	✓	✓	✓	✓	✓
Control III	✓	✓		✓	✓
Coverage™ Germicidal Surface Wipes	✓	✓	✓	✓	✓
Coverage™ Plus Germicidal Surface Wipes	✓	✓	✓	✓	
Coverage™ Spray			✓		
Dispatch Wipes	✓	✓		✓	✓
Dynacide PA	✓	✓		✓	✓
Echo Clean Lingettes		✓		✓	✓
End-Bac® II	✓	✓		✓	✓
Endosporine	✓	✓	✓	✓	✓
Envirocide®			✓		
Gluco Chlor Wipes	✓	✓		✓	✓
Hydrogen Peroxide (3%)	✓	✓	✓	✓	✓
Incidin Plus (1%)	✓	✓	✓	✓	✓
Incidin Plus (3%)	✓	✓	✓	✓	✓

Table 7: Disinfectant Compatibility with SonoSite Specialty Transducers (continued)

Disinfectant	C8 C8e C8x	ICT ICTx	ICTn	SLA SLAx	SLT
Korsolex® Extra	✓	✓		✓	✓
Madacide 1		✓	✓	✓	✓
Medzyme	✓	✓	✓	✓	✓
Metricide OPA Plus	✓	✓	✓	✓	✓
MetriSponge	✓	✓	✓	✓	✓
PerCept RTU Wipes™		✓	✓	✓	✓
Rely+On™ PeraSafe™					✓
Sani-Cloth® Bleach Wipes	✓	✓	✓	✓	✓
Sani-Cloth® HB	✓	✓	✓	✓	✓
Sani-Cloth® Plus	✓	✓		✓	✓
Sekusept®	✓	✓		✓	✓
Staphene®	✓	✓		✓	✓
Steranios (2%)	✓	✓	✓	✓	✓
Steranios (20%)					✓
T-Spray™		✓	✓	✓	✓
T-Spray II™	✓	✓	✓	✓	✓
Trigene Advance Wipes	✓	✓	✓	✓	✓
Tristel Duo	✓				
Tristel Solo	✓	✓	✓	✓	✓
Tristel Wipes					✓
Virex® 256			✓		

Table 7: Disinfectant Compatibility with SonoSite Specialty Transducers (continued)

Disinfectant	C8 C8e C8x	ICT ICTx	ICTn	SLA SLAx	SLT
Virex® TB			✓		
Virox® 5 Wipes	✓	✓	✓	✓	✓
Wet Task Wipes	✓	✓	✓	✓	✓
Wet Wipe Disinfection	✓	✓		✓	✓
Wex-Cide	✓	✓		✓	
Wipes Plus Disinfecting Wipes		✓	✓	✓	✓

✓ - Acceptable for use

Accessories

You can disinfect the exterior surface of the accessories using a recommended disinfectant.

Table 8: Disinfectant Compatibility with SonoSite Accessories

Disinfectant	Edge™ Stand	H-Universal™ Stand	V-Universal™ Stand	MDS Lite, MDS Lite II	Basic Stand, Footswitch, iLook® Stand, MDS®, MDSm, S Series™ Stand, ScanPack, SiteCharge® SiteStand®,	System battery	iLook® docking station	NanoMaxx® dock	PowerPack	PowerPark	ECG cable	SonoRemote™ control	Field Case (NanoMaxx® system)	Power Supply P06936, P08792, P09823, P09714	Triple Transducer Connect
Airkem A-33®	✓	✓	✓					✓	✓	✓			✓	✓	✓
Alcohol, Ethanol								✓							
Alcohol, Isopropyl (100%)	✓	✓	✓												✓
Alcohol, Isopropyl (70%)				✓	✓	✓	✓	✓			✓				
Anios Wipes	✓	✓	✓					✓	✓	✓		✓		✓	✓
Aniosyme DD1 (5%)	✓	✓	✓					✓	✓	✓				✓	✓
Anioxyde 1000	✓	✓						✓	✓	✓		✓		✓	✓
Asepti-HB	✓	✓						✓	✓	✓			✓	✓	✓
Asepti-Steryl®								✓	✓					✓	✓
Asepti-Wipe™ II	✓	✓						✓	✓	✓				✓	✓
Baxedin Pre-Op			✓					✓	✓						
CaviWipes™	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓
Cidalkan®								✓						✓	✓

Table 8: Disinfectant Compatibility with SonoSite Accessories (continued)

Disinfectant	Edge™ Stand	H-Universal™ Stand	V-Universal™ Stand	MDS Lite, MDS Lite II	Basic Stand, Footswitch, iLook® Stand, MDSe, MDSm, S Series™ Stand, ScanPack, SiteCharge® SiteStand®,	System battery	iLook® docking station	NanoMaxx® dock	PowerPack	PowerPark	ECG cable	SonoRemote™ control	Field Case (NanoMaxx® system)	Power Supply P06936, P08792, P09823, P09714	Triple Transducer Connect
Cidalkan® Wipes								✓							
Cidex®	✓	✓	✓					✓	✓	✓				✓	✓
Cidex® OPA	✓	✓	✓					✓	✓	✓			✓	✓	✓
Cidex® Plus	✓	✓	✓					✓	✓	✓				✓	✓
Clinell Universal Sanitising Wipes	✓	✓	✓					✓	✓					✓	✓
Clorox® Germicidal Wipe								✓	✓					✓	
Control III	✓	✓	✓					✓	✓	✓				✓	✓
Coverage™ Germicidal Surface Wipes	✓	✓	✓					✓						✓	✓
Coverage™ Plus Germicidal Surface Wipes								✓						✓	✓
Coverage™ Plus Wipes												✓			
Coverage™ Spray	✓	✓	✓					✓	✓					✓	✓
Coverage™ Wipes									✓			✓			
DisCide Ultra Disinfecting Towelettes								✓	✓					✓	✓

Table 8: Disinfectant Compatibility with SonoSite Accessories (continued)

Disinfectant	Edge™ Stand	H-Universal™ Stand	V-Universal™ Stand	MDS Lite, MDS Lite II	Basic Stand, Footswitch, iLook® Stand, MDSe, MDSm, S Series™ Stand, ScanPack, SiteCharge® SiteStand®,	System battery	iLook® docking station	NanoMaxx® dock	PowerPack	PowerPark	ECG cable	SonoRemote™ control	Field Case (NanoMaxx® system)	Power Supply P06936, P08792, P09823, P09714	Triple Transducer Connect
Dynacide PA	✓	✓	✓					✓	✓	✓				✓	✓
Echo Clean Lingettes								✓	✓	✓				✓	✓
End-Bac® II								✓	✓					✓	✓
Endosporine	✓	✓	✓					✓	✓	✓		✓		✓	✓
Envirocide®	✓	✓	✓					✓	✓	✓				✓	✓
Gluco Chlor Wipes								✓	✓					✓	
Hydrogen Peroxide (3%)	✓	✓	✓					✓	✓	✓		✓ ¹	✓	✓	✓
Incidin Plus (1%)	✓	✓	✓					✓	✓	✓				✓	✓
Incidin Plus (3%)	✓	✓	✓					✓	✓	✓				✓	
Madacide 1								✓	✓					✓	
Medzyme	✓	✓	✓					✓	✓	✓				✓	
Metricide OPA Plus	✓	✓	✓					✓	✓	✓		✓	✓	✓	✓
Mikrobac Forte 2%	✓														
Mikrobac Tissues	✓														
Nu-Cidex								✓		✓					✓

Table 8: Disinfectant Compatibility with SonoSite Accessories (continued)

Disinfectant	Edge™ Stand	H-Universal™ Stand	V-Universal™ Stand	MDS Lite, MDS Lite II	Basic Stand, Footswitch, iLook® Stand, MDSe, MDSm, S Series™ Stand, ScanPack, SiteCharge® SiteStand®,	System battery	iLook® docking station	NanoMaxx® dock	PowerPack	PowerPark	ECG cable	SonoRemote™ control	Field Case (NanoMaxx® system)	Power Supply P06936, P08792, P09823, P09714	Triple Transducer Connect
Oxivir TB Wipes								✓							
PerCept RTU Wipes™	✓	✓	✓					✓	✓	✓		✓	✓	✓	
Rely+On™ PeraSafe™	✓	✓						✓	✓	✓		✓		✓	✓
Sani-Cloth® AF	✓														✓
Sani-Cloth® Bleach Wipes								✓	✓					✓	
Sani-Cloth® HB								✓	✓			✓		✓	
Sani-Cloth® Plus	✓	✓		✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Simple Chlor								✓							
Sporicidin®	✓	✓	✓					✓	✓	✓		✓		✓	✓
Staphene								✓	✓					✓	
Steranios 2%	✓	✓	✓					✓	✓	✓				✓	✓
Steranios 20%	✓	✓	✓					✓	✓	✓				✓	✓
Super Sani-Cloth®					✓	✓	✓	✓	✓		✓	✓ ¹		✓	✓
T-Spray™		✓	✓						✓					✓	
T-Spray II™	✓	✓	✓					✓	✓	✓		✓	✓	✓	✓

Table 8: Disinfectant Compatibility with SonoSite Accessories (continued)

Disinfectant	Edge™ Stand	H-Universal™ Stand	V-Universal™ Stand	MDS Lite, MDS Lite II	Basic Stand, Footswitch, iLook® Stand, MDSe, MDSm, S Series™ Stand, ScanPack, SiteCharge® SiteStand®,	System battery	iLook® docking station	NanoMaxx® dock	PowerPack	PowerPark	ECG cable	SonoRemote™ control	Field Case (NanoMaxx® system)	Power Supply P06936, P08792, P09823, P09714	Triple Transducer Connect
Trigene Advance Wipes								✓						✓	
Tristel Solo	✓	✓						✓	✓	✓		✓	✓	✓	✓
Tristel Wipes								✓						✓	
Virex® 256		✓						✓	✓	✓			✓	✓	
Virex® II 256	✓														
Virex® TB	✓	✓	✓					✓	✓					✓	
Virox® 5 Wipes	✓	✓	✓					✓	✓	✓		✓		✓	✓
Wet Task Wipes	✓	✓	✓					✓	✓					✓	
Wex-Cide	✓	✓	✓					✓	✓	✓		✓		✓	✓
Wipes Plus Disinfecting Wipes	✓	✓	✓					✓	✓					✓	✓
Wood Wyant Ultra Wipes	✓	✓	✓					✓							✓

✓ - Acceptable for use

1 Not for use on the microphone clip. (The microphone clip is detachable from the microphone.)

Cleaning transducers

Clean the transducer prior to disinfecting to remove any particulate matter or bodily fluids.

To clean a transducer

- 1 Disconnect the transducer from the ultrasound system.
- 2 Remove any transducer sheath.
- 3 Clean the surface using a soft cloth lightly dampened in a mild soap or detergent cleaning solution.
Apply the solution to the cloth rather than the surface.
- 4 Rinse with water or wipe with water-dampened cloth; then wipe with a dry cloth.

Edge, H-Universal, iLook, M-OB/GYN Office, M-Turbo, MicroMaxx, NanoMaxx, S Series, SiteCharge, SiteStand, SonoRemote, SonoSite, SonoSite 180PLUS, TITAN, V-Universal, and the SonoSite logo are registered (in some jurisdictions) and unregistered trademarks owned by SonoSite, Inc.

The SonoSite ultrasound system(s) referenced in this document may be covered by one or more of the following U.S. patents: 5722412, 5817024, 5893363, 6135961, 6203498, 6364839, 6371918, 6383139, 6416475, 6447451, 6471651, 6569101, 6648826, 6575908, 6604630, 6817982, 6835177, 6962566, 7169108, 7449640, 7534211, 7549961, 7588541, 7591786, 7604596, 7643040, 7686766, 7694814, 7727153, 7740586, 7804970, 7809400, 7819807, 7841575, 7849250, 7867168, 7883276, D456509, D461895, D509900, D538432, D544962, D558351, D559390, D591423, D592750, D592760, D625014, D625015, and by the following counterpart foreign patents: AU727381, AU730822, CA2372152, CA2372158, CA2373065, CN ZL 97113678.5, CN ZL 98106133.8, CN ZL 98108973.9, CN ZL 200830007734.8, DE60021552.0, DE60029777.2, DE60034670.6, DE69730563.5, DE6980539.6, DE69831698.3, DE60 2004 23 816.3-08, FR0815793, FR0875203, FR0881492, FR1175713, FR1180970, FR1589878, GB0875203, GB0881492, GB1175713, GB1180970, GB1180971, GB1589878, IT0815793, IT0881492, IT1175713, IT1589878, KR528102, KR532359, NO326202, NO326814, NZ542968, RCD000897368-0001, SP0815793, SP0881492, SP1589878. Patents pending.

SonoSite, Inc.

21919 30th Drive SE
Bothell, WA 98021 USA
T: 1-888-482-9449 or
1-425-951-1200
F: 1-425-951-1201

SonoSite Ltd

Alexander House
40A Wilbury Way
Hitchin, Herts
SG4 0AP UK
T: +44-1462-444800
F: +44-1462-444801

P06703-08 11/2011
Copyright 2011 by
SonoSite, Inc.
All rights reserved.

